

Immunoassay

REF

CMB0201 / CMB0202 / CMB0203 / CMB0204 / CMB0205

50 tests*1 / 100 tests*1 / 100 tests*2 / 100 tests*5 / 50 tests*2

CEA CLIA Microparticles

This assay is based on a chemiluminescent microparticle immunoassay (CLIA Microparticles) for the quantitative determination of CEA (carcinoembryonic antigen) in human serum.

All trademarks are properties of their respective owners.

Key to Graphical Symbols Used

LOT

batch code



use by



manufacturer



contains sufficient for <n> tests

IVD
in vitro diagnostic medical device


temperature limitation

REF

catalogue number



consult instructions for use

EC REP

authorized representative in the European Community

EC REP

 OBELIS S.A.
 Bd. Général Wahis, 53
 1030 Brussels
 Belgium

 AUTOBIO DIAGNOSTICS CO., LTD.
 No.87 Jingbei Yi Road
 National Eco & Tech Development Area
 Zhengzhou
 China
 450016


For any technical assistance please contact us in English at:
 Email: customerservice@autobio.com.cn
 Contact your local dealers for all product-related questions in your local language

Introduction

CEA (carcinoembryonic antigen) is a cell-surface 200kD glycoprotein. It was first described by Gold and Freedman in 1965 as a complex immunoreactive glycoprotein found in epithelial adenocarcinomas of the colon and fetal colon.^{1,2} Increased levels of CEA are observed in more than 30% of patients with cancer of the lung, liver, pancreas, breast, colon, head or neck, bladder, cervix, and prostate. Elevated plasma levels are related to the stage and extent of the disease, the degree of differentiation of the tumor, and the site of metastasis. Its main use is in the monitoring of cancer patients after surgery, chemotherapy or radiotherapy.³ Successful removal of the tumor is usually followed by a decrease in the concentration of circulating CEA,⁴ whereas recurrence of the primary tumor or metastasis is accompanied by increasing concentrations.⁵ Elevated serum levels of CEA may be found in a variety of benign and malignant conditions other than carcinoma of the colon. Conditions that may cause elevated levels of CEA include hepatic cirrhosis, hepatitis, heavy smoking, bronchitis, pancreatitis, gastritis, inflammatory bowel disease and renal disease.^{6,7}

Measurement Principle

This assay is based upon the one-step sandwich method. The sample, CEA antibody coated microparticles and enzyme labeled Anti-CEA are combined. CEA present in the sample is allowed to react simultaneously with the two antibodies, resulting in the CEA being sandwiched between the microparticles and enzyme-linked antibodies. After washing, a complex is generated among the microparticles, the CEA within the sample and enzyme-linked antibodies by immunological reactions. Chemiluminescent Substrate is added and the complex catalyzes substrate, resulting in a chemiluminescent reaction. The resulting chemiluminescent reaction is measured as RLU. The RLU is proportional to the amount of CEA in the samples.

Materials provided


1. Calibrators

6 vials each containing 1.0 mL of Calibrator A through F. The matrix is PBS (phosphate buffered saline) buffer containing casein. Contains a selection of preservatives.

Calibrators provided ready to use.

2. Reagent pack

Reagent pack provided ready to use.

	50*1	100*1	100*2	100*5	50*2
Microparticles Solution	1.2mL*1	2.3mL*1	2.3mL*2	2.3mL*5	1.2mL*2
Enzyme Conjugate	5.5mL*1	11.0mL*1	11.0mL*2	11.0mL*5	5.5mL*2

● Microparticles Solution

CEA antibody coated microparticles in PBS buffer containing BSA (bovine serum albumin). Contains a selection of preservatives.

● Enzyme Conjugate

Horse radish peroxidase labeled mouse monoclonal Anti-CEA in Tris buffer containing bovine serum. Contains a selection of preservatives.

Assay Analyzers on which the kit can be used

- AutoLumo A2000 Plus
- AutoLumo A2000 Plus B
- AutoLumo A1000

The chemiluminescent microparticle immunoassay (CLIA Microparticles) is intended for use on Assay Analyzer which is AutoLumo A2000 Plus,

CEA CLIA Microparticles

AutoLumo A2000 Plus B or Autolumo A1000.

Materials Required but not Provided

1. Assay Analyzer
2. Reaction vessel(s) for sample and reagent reaction
3. Sample tube(s) or cup(s) for sample containing
4. Diluent Universal
5. Chemiluminescent Substrate
6. System Wash for washing the pipetting needles
7. Wash Buffer used in the washing procedure
8. Distilled or deionized water

Metrological Traceability of Calibrators

The measurand or analyte in the CEA calibrators is traceable to the material purchased from NIBSC (National Institute for Biological Standards and Control), the NIBSC code: 73/601.

Warnings and Precautions

1. For professional use only.
2. Follow the instruction for use carefully. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this instruction for use.
3. Refer to the material safety data sheet and product labeling for any chemical hazards that may be present in this assay.
4. Handle the potentially contaminated materials and wastes safely according to local requirement.
5. This assay contains materials of animal origin. Bovine components originate from countries where bovine spongiform encephalopathy (BSE) has not been reported.
6. Do not smoke, drink, eat or use cosmetics in the working area.
7. Wear protective clothing and disposable gloves when dealing with samples and reagents. Wash hands after operations.
8. Conduct the assay away from bad ambient conditions. e. g. ambient air containing high concentration corrosive gas such as sodium hypochlorite acid, alkaline, acetaldehyde and so on, or containing dust.
9. Do not use reagents beyond the labeled expiry date.
10. Do not mix or use components from kits with different batch codes.
11. When storing the calibrators, be certain the vials are securely sealed.
12. Ensure the microparticles are resuspended before loading it on the analyzer.
13. Avoid foam formation in all reagents and sample types (samples, calibrators and controls).
14. Do not substitute any reagent in this kit from other manufacturers or other lots.
15. When any damage to the protective packaging or any change of analytical performance is observed, do not use the kit.

Storage

1. Store the kit at 2-8°C. Do not freeze. Avoid strong light. When stored as directed, all reagents are stable until the expiration date.
2. Refrigerate the reagent pack at 2-10°C for a minimum of 2 hours prior to use.
3. Store the unsealed reagents pack upright on the analyzer or 2-10°C for a maximum of 28 days. After 28 days, the reagent pack must be discarded. Once they are removed from the analyzer, store them at 2-10°C in an upright position.
4. Seal and return the remaining calibrators at 2-8 °C immediately after the experiment, under which conditions the stability will be retained for 1 month, for longer use, store opened calibrators in aliquots and freeze at -20°C. Avoid multiple freeze-thaw cycles, do not freeze-thaw

more than 3 cycles.

Sample

1. Collect serum samples in accordance with correct medical practices.
2. Do not use heat-inactivated samples. Do not use sodium azide preservative in samples.
3. Do not use samples with obvious microbial contamination.
4. Sediments and suspended solids in samples may interfere with the test result which should be removed by centrifugation. Ensure that complete clot formation in serum samples has taken place prior to centrifugation. Some samples, especially those from patients receiving anticoagulant or thrombolytic therapy, may exhibit increased clotting time. If the sample is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results. Be sure that the samples are not decayed prior to use.
5. Prior to shipment, it is recommended that samples be removed from the clot, serum separator or red blood cells.
6. Use caution when handling patient samples to prevent cross contamination. Use of disposable pipettes or pipette tips is recommended.
7. Insufficient processing of sample or disruption of the sample during transportation may cause depressed results.
8. Avoid grossly hemolytic, lipemic or turbid samples.
9. Cap and store the samples at 18-25°C for no more than 8 hours, for longer use samples should be capped and stored at 2-8°C up to 48 hours. Or freeze the samples that need to be stored or transported for more than 48 hours at -20°C. Avoid multiple freeze-thaw cycles. Mix thawed samples thoroughly by low speed vortex or by inverting 10 times. Visually inspect the samples, if layering or stratification is observed, continue mixing until samples are visibly homogeneous. After thawing, bring to room temperature and mix well by gently shaking.
10. Centrifuge the thawed samples containing red blood cells or particulate matter, or which are hazy or cloudy in appearance prior to use to ensure consistency in the results.
11. Note that interfering levels of fibrin may be present in samples that do not have obvious or visible particulate matter.
12. If proper sample collection and preparation cannot be verified, or if samples have been disrupted due to transportation or sample handling, an additional centrifugation step is recommended. Centrifugation conditions should be sufficient to remove particulate matter.
13. For optimal results, inspect all samples for bubbles. Remove bubbles with a tip prior to analysis. Use a new tip for each sample to prevent cross contamination.

Measurement Procedure

1. **Check the consumable materials**
 - Verify adequate volume of consumable materials is present prior to running the test.
 - Refer to the Assay Analyzer's operation manual.
2. **Load the kit**
 - Mix contents of new (unpunctured) reagent packs by gently inverting pack several times before loading on the analyzer. Avoid foam formation in all reagents. Don't invert the open (punctured) packs. If necessary, shake gently to mix horizontally after the first loading.
 - Read the bar code on the reagent pack automatically to obtain the required parameters for the test.
 - If the bar code cannot be read in exceptional cases, they can be recognized manually.
 - Refer to the Assay Analyzer's operation manual.
3. **Order tests**

- Place the sample tubes or cups on the sample rack, 25 μ L of samples for each test. But consider the sample container and 150 μ L of system dead volumes, which can be refer to the appropriate Assay Analyzer manuals for the minimum sample volume required.
- Load the sample rack and input the sample information on the system software interface.
- Select "run" to start the test, the analyzer automatically operates tests. It performs the following functions:
 - Moves the sample to the set point
 - Loads a reaction vessel into the process path
 - Aspirates and transfers sample into the reaction vessel
 - Adds Microparticles Solution and Enzyme Conjugate to the reaction vessel
 - Mixes, incubates and washes the reaction mixture
 - Adds Chemiluminescent Substrate
 - Measures chemiluminescent emission to determine the quantity of CEA in the sample
 - Discards the used reaction vessel
 - Calculates the result
- Refer to the Assay Analyzer's operation manual.

4. Calibrate the curve

- Analyzer can read the bar code on the reagent pack automatically to obtain the required parameters for the test.
- If the bar code cannot be read in exceptional cases, they can be recognized manually.
- Transfer the calibrators into the sample tubes or cups and place them on the sample rack. Conduct duplicate detection on the system.
- Load the sample rack and input calibrators' information on the system software interface.
- Select "run" to start the test and generate the calibration curve, calibration is required every 28 days.
- Once a calibration curve is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - Controls are out of range after repeated measurements
 - A reagent kit and Chemiluminescent Substrate with new batch code is used
 - Beyond the expiration date of a calibration curve
 - Important parts of the analyzer are replaced or repaired
- Refer to the Assay Analyzer's operation manual.

5. Dilute the sample

Samples with a CEA exceeding 1000 ng/mL may be diluted via the program of the analyzer. Diluent Universal is used to dilute the samples. The software takes the dilution into account when reporting the result.

- The concentration of the sample after dilution should not be less than 50 ng/mL.

Measurement Results

The sample test results are determined automatically by the system software. The amount of CEA in the samples is determined from the measured light production by means of the stored calibration data. Refer to the Assay Analyzer's operation manual on reviewing sample results.

Control Procedure

The recommended control requirement for this assay is to purchase control materials separately and test them together with the samples within the same run. The result is valid if the control values fall within the acceptable ranges. When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results may be invalid and may require retesting. Assay recalibration

may be necessary. It is recommended that each laboratory establish its accepted range to ensure proper test performance.

Limitations of the Procedure

1. This assay is intended as an aid for the clinical diagnosis. Conduct this assay in conjunction with clinical examination, patient's medical history and other test results.
2. If the results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
3. Heterophilic antibodies and rheumatoid factors in samples may interfere with test results. Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed. Additional information may be required for diagnosis. This kind of samples is not suitable to be tested by this assay.
4. Performance of this test has not been established with neonatal samples.
5. CEA level in the same range as healthy individuals. Elevations in circulating CEA levels may be observed in smokers as well as in patients with nonmalignant disease. For these reasons, a serum CEA level, regardless of value, should not be interpreted as absolute evidence for the presence or absence of malignant disease. The CEA level should be used in conjunction with information available from clinical evaluation and other diagnostic procedures. This CEA assay should not be used as a cancer screening test.
6. Patients who have received mouse monoclonal antibodies for either diagnosis or therapy can develop HAMA (human Anti-mouse antibodies). HAMA can produce either falsely high or falsely low values in immunoassays which use mouse monoclonal antibodies. Additional information may be required for diagnosis.
7. This test measures concentrations within the range of 0.5-1000 ng/mL. If CEA concentrations above the measuring range to be expected, it is recommended to dilute samples with Diluent Universal, the recommended dilution is 1:3 of this test, under this condition, allowing samples to be up to approximately 4000 ng/mL.

Biological Reference Interval

A normal value of <5 ng/mL (95th percentile) was obtained by testing serum samples from 704 individuals defined as normal by clinician. It is recommended that each laboratory establish its own normal range which may be unique to the population it serves depending upon geographical, patient, dietary, or environmental factors.

Performance Characteristics

1. Measurement Precision

3 samples were assayed in duplicate, twice per day across 20 testing days. Data from this study are summarized in the following table.

Sample	n	Mean (ng/mL)	Within-run	Total
			%CV	%CV
1	80	5.81	3.45	4.97
2	80	10.87	1.88	4.94
3	80	98.84	2.09	4.84

*Representative data; results in individual laboratories may vary from these data.

2. Analytical Sensitivity

Limit of Blank: 0.5 ng/mL

Limit of Detection: 1.5 ng/mL

Limit of Quantitation: 1.8 ng/mL with a coefficient of variation of ≤ 20 %

3. Analytical Specificity

Cross reaction: the following substances and concentrations were tested and found no cross reaction with the test.

Substances	Concentration
AFP	500 ng/mL
CA19-9	500 U/mL
CA125	400 U/mL
CA15-3	500 U/mL

Interference: No interference with 30 mg/dL of bilirubin, 400 mg/dL of hemoglobin, 3000 mg/dL of triglycerides.

4. Measurement Accuracy by Correlation

A comparison study was performed where samples were tested using this assay and a CEA reference assay. Data were analyzed and are summarized in the following table.

Correlation Method	Number of Samples	Intercept	Slope	Correlation Coefficient
Linear Regression	457	1.3083	0.9734	0.982

5. High Dose Hook Effect

A sample spiked with CEA up to 100,000 ng/mL was determined, the concentration result obtained was ≥ 1000 ng/mL.

Literature References

1. GOLD P, FREEDMAN SO. DEMONSTRATION OF TUMOR-SPECIFIC ANTIGENS IN HUMAN COLONIC CARCINOMATA BY IMMUNOLOGICAL TOLERANCE AND ABSORPTION TECHNIQUES. J. Exp. Med. 1965;121:439-462.
2. Gold P, Freedman SO. Specific carcinoembryonic antigens of the human digestive system. J. Exp. Med. 1965;122(3):467-481.
3. Smith AM, Macdonald DJ. Enzyme immunoassay for carcinoembryonic antigen. Clin. Chem. 1983;29(12):2019-2021.
4. Lokich JJ, Zamcheck N, Lowenstein MW. Sequential carcinoembryonic antigen levels in the therapy of metastatic breast cancer: a predictor and monitor of response and relapse. Ann. Intern. Med. 1978;89(6):902-906.
5. Falkson HC, Falkson G, Portugal MA, van der Watt JJ, Schoeman HS. Carcinoembryonic antigen as a marker in patients with breast cancer receiving postsurgical adjuvant chemotherapy. Cancer. 1982;49(9):1859-1865.
6. Torosian MH. The clinical usefulness and limitations of tumor markers. Surg Gynecol Obstet. 1988;166(6):567-579.
7. Sugarbaker PH. Role of carcinoembryonic antigen assay in the management of cancer. Adv Immun Cancer Ther. 1985;1:167-193.